

PSU NMR Facility


Acquiring an ^1H NMR spectrum


(Quick instructions for Topspin 1.3, Bruker DRX-400)

If you need assistance, please contact:

- 1) Dr. Tapas Mal -- Room #8 - tkm9@psu.edu; tapas1269@gmail.com;
Ph# office: (814) 865-0941; cell: (617) 803-2411
- 2) Dr. Carlos Pacheco -- Room#14 - cnp14@psu.edu; cpacheco223@gmail.com;
Ph# office: (814)863-1182; cell: (609) 240-5957
- 3) Dr. Debashish Sahu -- Room#10 - desahu@psu.edu; dzs12@psu.edu;
debashish.sahu@gmail.com

Place your sample in the spinner and carefully adjust the NMR tube depth in the gauge.


- 1) **Turn on** the lift air using BSMS (or type **ej** <enter> at the command line).
 - a. If BSMS window is missing, type **bsmsdisp** <enter> at the command line.
 - b. **Always wait for the sound of the lift air at full power before placing your sample on top of the magnet!!!**
- 2) **Turn off** lift air from BSMS (or type **ij** <enter> at the command line).
- 3) Create a new data set by typing **new** <enter> at the command line or click on  at the toolbar (or type **edc** <enter> at the command line)



NAME -- (filename for your compound)
 EXPNO -- suggestion: Set 1 for a ^1H , 2 for a ^{13}C
 PROCNO -- suggestion: set to 1
 DIR -- set to /opt/topspin/
 USER -- should be <Your_Professor last name>
 SOLVENT -- your deuterated solvent from the pull-down menu
 EXPERIMENT -- "Use current params"
 TITLE -- the title of your sample (can be as descriptive as needed)

Read a good shimfile with **rsh**<enter>: choose the **stdshims** file from the list, which is the standard shimfile name for all NMR spectrometers.

- 4) Spin the sample at 20 Hz: type **ro** <enter> at the command line and click on **Start rotation** at the window:



- 5) Type **lock** <enter> at the command line and choose the appropriate solvent from the lock table.
- 6) If the lock window is missing, type **lockdisp** <enter> at the command line
- 7) After the lock process is finished, **shim z1 and z2**.
- 8) Type **rpar PROTON all** <enter> at the command line.
- 9) Type **getprosol** <enter> at the command line.
- a. For instruments running **Xwinnmr**, type **eda** <enter> and set prosol to **TRUE**.
- 10) Type **rga** <enter> at the command line
- a. Adjust the number of scans - **NS** - as appropriate if necessary; default for ^1H is 16, for ^{13}C is 1024).
- 11) Type **zg** <enter> at the command line to start the experiment.
- 12) For instruments running **Xwinnmr**, type **acqu** <enter> at the command line if you do not see your FID (raw NMR data).
- 13) Your NMR data is automatically saved **at the end** of the experiment.

DATA PROCESSING

It is highly recommended to carry out NMR data processing on the *off-line* PC at the NMR Facility, room #10, where you can also back your NMR data up

****Bruker is now offering Topspin software free for Academia: the software can be downloaded at:**

<https://www.bruker.com/service/support-upgrades/software-downloads/nmr/free-topspin-processing/free-topspin-download.html>

Optionally, you may use the MNova software (the NMR Facility has 50 licenses) – consult the NMR staff to obtain your license and instructions for downloading and activation

- 14) Type **ef; apk** <enter> at the command line to perform Apodization (**e**), Fourier (**f**) Transformation, and automatic phase adjustment (**apk**).
- 15) If needed, further optimize the phase of the spectrum manually, **using the phase correction icon**.


Click on the phase icon:




Left-click **0**, hold the mouse button to correct phase **at** the vertical cursor.

Left-click **1**, hold the mouse button to correct phase **away** from the cursor.

Click on the SAVE icon to exit the phasing window.




- 16) Calibrate your spectrum using a known chemical shift (residual solvent peak or TMS). Click on the icon  at the toolbar. Set the reference and return.

- 17) Type **abs** to flatten the baseline. There is a manual baseline flattening:




18) Integrate the spectrum manually with  -- This takes you into a subroutine:

Integrate the spectrum, save the integrals before exit the window.



19) Peak picking can be done with the icon  and drawing a green box around the peaks

Save the peaks before exiting the window (using the same icon as the phasing and integration).



20) Type **plot** <enter> to open the plot editor (for instruments running TopSpin) and then print.

- a. For instruments running Xwinnmr type **view** <enter> to print preview and type **plot** <enter>.

21) Print.

22) Save your data on the PC at Room #10

- a. Save your NMR data **ALSO** on a thumb drive (recommended).
 - i. Email the NMR data (as .zip file) is also possible at the off-line PC.